# Lysis buffer for DirectAce 組織溶解液 (DirectAce 用)

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# Lysis buffer for DirectAce

## I. Description

Lysis buffer for DirectAce is a DNA extraction buffer optimized as a processing reagent for DirectAce qPCR Mix plus ROX Tube (Code No. 318-07751), Real time PCR reagent for detection of fluorescence-labeled probe.

## II. Storage Temperature

Room temperature

## III. Contents

Lysis buffer for DirectAce 100 ml X 1 本

## IV. Storage Conditions

20 mM Tris-HCl (pH 8.0)

400 mM NaCl

5 mM EDTA

0.3% SDS

## V. Cautions

- It is possible to change the amount ratio of sample and Lysis buffer depending on a kind and a condition of the sample.
- If necessary, centrifuge the sample after incubation.

## VI. Example of use

- ① using Food mill  
Crush 5 g of sample with 20 ml of Lysis buffer for DirectAce by using food mill.  
Incubate the sample for 20 minutes at RT.  
Dilute the supernatant to 2 times or more by water or Lysis buffer for DirectAce.
- ② extract from ground sample  
Add 400  $\mu$ l of Lysis buffer for DirectAce to 50 mg of sample, and vortex it.  
Incubate the sample for 20 minutes at RT.  
Dilute the supernatant to 2 times or more by water or Lysis buffer for DirectAce.

### ③ using pestle

Crush 50 mg of sample with 400  $\mu$ l of Lysis buffer for DirectAce by using pestle, and vortex it.

Incubate the sample for 20 minutes at RT.  
Dilute the supernatant to 2 times or more by water or Lysis buffer for DirectAce.

### <References>

Mano, J.; Hatano, S.; Futo, S.; Minegishi, Y.; Ninomiya, K.; Nakamura, K.; Kondo, K.; Teshima, R.; Takabatake, R.; Kitta, K.; Development of direct real-time PCR system applicable to a wide range of foods and agricultural products., Food Hyg. Saf. Sci. 2014, 55, 25-33.

For research use only.

Not for use in diagnostic or therapeutic purpose.