Uracil-DNA Glycosylase (UNG)

I. Description

E. coli Uracil-DNA Glycosylase (UNG) removes uracil from DNA by hydrolysis of the N-glycosidic bond between uracil base and deoxyribose sugar, leaving an apyrimidinic site in the DNA.

UNG prevents carry-over of DNA in PCR products incorporating dU-nucleotides. To eliminate carryover of PCR products from previous reactions, use a qPCR master-mix reagent that contains dUTP, and add the UNG treatment step prior to starting PCR.

II. Storage

Store at -20°C

III. Source

Escherichia coli

IV. Components

Product	Conc.	100 ng
Uracil-DNA Glycosylase (UNG)	1 ng / μL	100 μL x 1

V. Purity

- 1 ng of the enzyme and 1 μg of substrate DNA (normal DNA not containing uracil) were reacted at 37°C for 16 hr in 50 μL of the reaction volume and then agarose gel electrophoresis was performed. The result indicated no degradation of DNA.
- 1 ng of the enzyme and 1 μg of substrate RNA were reacted at 37°C for 2 hr in 50 μL of the reaction volume and then agarose gel electrophoresis was performed. The result indicated no degradation of RNA.

VI. Notes

- Too much UNG may inhibit the PCR. Suitable amount of the enzyme is commonly in the range between 0.1-1 μL in a 50 μL reaction volume.
- Do not use this product for One-step RT-PCR application as the UNG would active on uracil in cDNA among 55°C .

VII. Usage

To prevent carryover contamination in your qPCR, use a qPCR master-mix reagent that contains dUTP, and add UNG in reaction mix and incubate at 40°C for 2 min prior to starting qPCR reactions.

An example of using GeneAce Probe qPCR Mix II (Code No. 313-08823):

Component	Volume	Volume
UNG (1 ng/ μL)	0.5 μL	0.25 μL
2 x GeneAce Probe qPCR Mix II	25.0 μL	12.5 μL
25 μM Forward Primer	1.0 μL	0.5 μL
25 μM Reverse Primer	1.0 μL	0.5 μL
10 μM TaqMan [®] Probe ^{*1)}	1.0 μL	0.5 μL
Template DNA	5.0 μL	2.5 μL
d.d.H ₂ O	up to	50.0 μL 25.0 μL

UNG and PCR conditions

40°C	2 min.) 45 cycles
95°C	10 min.	
95°C	30 sec.	
60°C	1 min.	

UNG treatment
Heat-inactivate UNG ^{*2)}

*1) TaqMan[®] is registered trademarks of Roche Molecular Systems, Inc.

*2) When GeneAce Probe qPCR Mix II is used, the PCR must start with initial 10 minutes incubation at 95°C to activate the chemically modified hot-start Taq DNA polymerase.

VIII. Related product

A set of qPCR reagent and UNG

Code No.	Product	Size
314-09051	GeneAce Probe qPCR Mix II with UNG	1 Set (200 rxns)

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