



# GM Maize Detection MIR162 Oligonucleotide Set

---

**Code No. 315-08361** (1 Set)

**Storage:**

Store at -20°C and protect from light.

**Component:**

Primer pair (25 μM each) 180 μL

Probe (10 μM) 180 μL

**Detection of:**

MIR162

**Reference:**

1. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.
2. Takabatake R, Masubuchi T, Futo S, Minegishi Y, Noguchi A, Kondo K, Teshima R, Kurashima T, Mano J, Kitta K.: *Food Hygiene and Safety Science (Shokuhin Eiseigaku Zasshi)*, **55** (5), 205 (2014)

**Probe & primer name:**

MIR162-1 [ MIR162-f1 & MIR162-r1 ]

MIR162-p1 (FAM/TAMRA)



# GMトウモロコシ系統別 DNA MIR162 オリゴヌクレオチドセット

Code No. 315-08361 (1 Set)

**保存:**

Store at -20°C(遮光)

**内容:**

Primer pair (25 μM each) 180 μL

Probe (10 μM) 180 μL

**検出:**

MIR162

**プローブ&プライマー名:**

MIR162-1 [ MIR162-f1 & MIR162-r1 ]

MIR162-p1 (FAM/TAMRA)

**参考文献:**

1. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」
2. Takabatake R, Masubuchi T, Futo S, Minegishi Y, Noguchi A, Kondo K, Teshima R, Kurashima T, Mano J, Kitta K.: *Food Hygiene and Safety Science (Shokuhin Eiseigaku Zasshi)*, **55** (5), 205 (2014)

本品は、試薬(試験研究用)として販売しているものです。  
医薬品の用途には使用しないでください。