



GMO Detection No Template Control -CoIE1/TE-

Code No. 318-06031 (90 μ L x 2)

Storage:

Store at -20°C

Component:

CoIE1 DNA/TE Solution 90 μ L x 2
(CoIE1 DNA 5 ng/ μ L)

Usage example:

Use 2.5 μ L of this product in each 25 μ L reaction as non-template control (e.g. Add 2.5 μ L of this product to 22.5 μ L of PCR mixture.)

Reference:

1. "JAS Analytical Handbook, Manual of Inspection and Analysis for GM-Foods"; in Japanese, by IAI Center for Food Quality, Labeling, and Consumer Services of Japan.
2. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.



GMO 共通 No Template Control –ColE1/TE–

Code No. 318-06031 (90 μ L x 2)

保存:

Store at -20°C

内容:

ColE1 DNA/TE 溶液 90 μ L x 2
(ColE1 DNA 5 ng/ μ L)

使用例:

PCR 系 25 μ L の場合、ブランク試料液として本品を 2.5 μ L 使用する (e.g. PCR 混合液 22.5 μ L に本品を 2.5 μ L 添加する。)

参考文献:

1. 農林水産消費技術センター「JAS 分析試験ハンドブック 遺伝子組換え食品検査・分析マニュアル」
2. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」

本品は、試薬 (試験研究用) として販売しているものです。
医薬品の用途には使用しないでください。