



GM Soybean (LLS) Detection LLS Oligonucleotide Set

Code No. 312-07391 (1 Set)

Storage:

Store at -20°C and protect from light.

Component:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 180 µL

Detection of:

LLS

Probe & primer name:

KVM175 & SMO001

TM031-Taq (FAM/TAMRA)

Reference:

1. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.
2. Takabatake R, Akiyama H, Sakata K, Onishi M, Koiwa T, Futo S, Minegishi Y, Teshima R, Mano J, Furui S, Kitta K.: *Shokuhin Eiseigaku Zasshi*, **52** (2), 100 (2011)

GM ダイズ(LLS)系統別 DNA LLS オリゴヌクレオチドセット

Code No. 312-07391 (1 Set)

保存:

Store at -20°C(遮光)

内容:

Primer pair (25 μM each) 180 μL

Probe (10 μM) 180 μL

検出:

LLS

プローブ&プライマー名:

KVM175 & SMO001

TM031-Taq (FAM/TAMRA)

参考文献:

1. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」
2. Takabatake R, Akiyama H, Sakata K, Onishi M, Koiwa T, Futo S, Minegishi Y, Teshima R, Mano J, Furuji S, Kitta K.: *Shokuhin Eiseigaku Zasshi*, **52** (2), 100 (2011)

本品は、試薬(試験研究用)として販売しているものです。
医薬品の用途には使用しないでください。